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LEONE MONOIMPLANTS FOR OVERDENTURE O-RING - INSTRUCTIONS FOR USE

Warnings and restrictions. The herewith furnished information is NOT intended for the patient, but exclusively for the medical professional, who will have to evaluate the contents in reference to every specific case to be treated, in order to furnish a suitable information to the patient about the benefits and eventual side effects or complications. The Leone monoimplants for overdenture O-ring must only be used by Dentists or Dental Surgeons experienced in implant techniques as well as in the related diagnostic methodologies, pre-operative planning and manufacture of prosthesis on implants. The hereunder instructions for use are finalized to underline the risks, the contraindications and the collateral effects of general character (a reference to all the possible circumstances of use is not possible) and well known at the moment of the placing on the market of the device. Nevertheless the Leone monoimplants for overdenture O-ring require a suitable knowledge of the Guidelines *Surgical Procedure and Prosthetic Procedure* available either on Leone Implantology Catalogue or on the web site www.leone.it. Considering that errors and deficiencies in patient evaluation, preoperative diagnosis, in the planning and the carrying out of the treatment may lead to implant loss, before use it is strongly recommended to attend specific training courses with the purpose to reach an high degree of propaedeutical knowledge and practice in the use. It is not ascribed to Leone the duty of monitoring the procedures of use of the products for the implant-prosthetic treatment; such procedures must integrally be reverted under the responsibility of the clinician. Leone takes no responsibility for the application of the device and its workmanship nor for possible surgical or prosthetic misuse of the device even neither for failure, adverse reactions or damages occurred to the patient or to the dentist following the application of the product. In any case the use of exclusively original components, accessories and instruments in conformity to the instructions is strongly recommended. **U.S.A. warning.** "Caution: federal law restricts this device to sale by or on the order of a physician." **Description of the product.** The Leone monoimplants for overdenture O-ring are medical devices made of titanium for surgical use (ISO 5832/2-3). The monoimplants for overdenture O-ring are provided with specific accessories which include components and instruments for surgery and prosthetic dentistry. The accessories and the instruments can be made of titanium, stainless steel and biocompatible plastic material. The monoimplants for overdenture O-ring are available in different sizes in order to match individual anatomical variability. **Intended use.** The Leone monoimplants for overdenture O-ring are indicated as therapy for the stabilization of removable prosthesis on lower arch. The Leone monoimplants for overdenture O-ring are designed to be surgically inserted in the bone structure of the mouth, exclusively in the mandible at the level of the area between the two foramina, as an anchorage system for total removable prosthesis. To properly support a removable prosthesis, four monoimplants have to be inserted. **Contraindications.** All the pathological conditions, both local and systemic, which contraindicate an intervention of oral surgery, have to be considered contraindications valid for the application of dental implants as well. As information only and provided that a specific clinical evaluation of the patient is strictly recommended, Leone monoimplants for overdenture O-ring should not be used in patients with recent myocardial infarct, haematological disorders, metabolic bone disorders, neoplastic pathologies diagnosed during the last 5 years, immunodepressed patients, previous bone radiotherapy or chemotherapy, patients with cardiovascular or metabolic problems which are not compensated, dismorphic patients, severe mental illness, xerostomia, diabetes or hyperglycaemia, transplanted patients, renal insufficiency, fibrous dysplasia, pharmacological treatments with anticoagulants, anticonvulsivants, immunosuppressants, steroids, acute or chronic infectious illnesses, pregnancy, nursing, insufficient healthy bone volume and impossibility or unavailability of the patient for the augmentation technique, skeletal growth still in action, scarce oral hygiene, parodontitis, candidosis and pathologies of the oral mucosa, allergic reaction to titanium, insufficient arch space, occlusion and/or articulation defects, para-functional habits (clenching, bruxism, etc.), poor patient's compliance and motivation, abuse of drugs, alcohol and tobacco, acute, sub-acute or chronic maxillary osteitis, radicular residuals. It is understood that the final choice to proceed or not with the treatment is an exclusive competence of the sanitary clinician. Do not insert the monoimplant in case of insufficient bone volume: at least 1.5 mm around the implant and 2 mm to the apical is recommended. The distance between one implant site and the adjacent one must necessarily be of at least 6 mm, to allow enough space for the positioning of the metal housings. The head of the monoimplant has to protrude from the gingiva for at least 1 mm to avoid a subsequent compression of the micro housings on the soft tissues. The monoimplants for overdenture O-ring should not touch, nor compress, or cut off the mandibular nerve. **Side effects.** Apart from the features of the utilized implant, the surgical application can be followed by edema, pain, haematomas, inflammations, limitation or alteration of the sensibility and oral functions, allergic reactions. Generally such effects are temporary and the patient has to receive a complete information about these side effects and about the necessity to give immediate communication to the physician when they would occur.

Complications and interactions. Referable complications related to the Leone monoimplants for overdenture O-ring are not known neither predictable. The surgical procedure, particularly if defectively performed, may cause hemorrhagic diseases, permanent paresthesia and/or dysesthesia, bone re-absorption at mandibular crest level, infectious processes, implant or prosthesis breakage, bone fractures, mucous ulcers, gingival hyperplasia, aesthetical problems, nervous lesions, inhalation or ingestion of parts or tools of the implant system. **Packaging and sterility.** The Leone monoimplants for overdenture O-ring are supplied gamma ray sterilized in glass containers. The product remains sterile if packaging is intact. When stored correctly, the sterility of the monoimplant is guaranteed up to the expiration date. Do not use the monoimplants after the expiration date printed on the label. Store the monoimplants in the original packaging, in a dry place, at room temperature and away from sun light. Leone takes no responsibility in case of re-sterilized monoimplants, regardless of the method of sterilization. Never use non-sterile or pre-used monoimplants. To open the monoimplant container and handle sterile products, please refer to explanatory illustrations. The accessories and the non-sterile single-use prosthetic components must be sterilized according to a validated procedure which is constantly monitored and follows ISO 17665-1 "Sterilization of health care products - moist heat - Part 1: requirements for the development, validation and routine control of a sterilization process for medical devices". It is recommended to sterilize the components with steam autoclave using the following parameters: temperature: 121°C (250°F), pressure: 1 atm, minimum exposure time: 20 minutes. The sterilization packaging must be appropriate for the device and conforming to EN 868-5 "Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods".

Identification and traceability. The monoimplants for overdenture O-ring are identified by means of the code and the lot number. This information is essential to assure the traceability of the product. The package is provided with an adhesive label showing these data and to be applied on the patient's record sheet. **Final precautions.** The patient should be informed about the precautions for the period after installation of the monoimplants, during and at the completion of treatment, in order to prevent complications, and variations in the efficiency of the device. For instance, heavy physical activity and occlusal trauma should be avoided. A good level of oral hygiene should be achieved and maintained. Periodical check-ups should be performed.

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