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
Figura 3: Distortion caused by stainless steel brackets detected by the 2 MRI sequences:

- 1, plastic trays, sagittal T1-weighted image;
- 2, stainless steel brackets, sagittal T1-weighted image;
- 3, plastic trays, axial GRE;
- 4, stainless steel brackets, axial GRE. [Elison, 2009]

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- Abbaszadeh K, Heffez LB, Mafee MF. Effect of interference of metallic objects on interpretation of T1-weighted magnetic resonance images in the maxillofacial region. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 2000;89:759–765.
- Arash Poorsattar-Bejeh Mira,b,* and Manouchehr Rahmati-Kamelc, Should the orthodontic brackets always be removed prior to magnetic resonance imaging (MRI)?, *J Oral Biol Craniofac Res.* 2016 May-Aug; 6(2): 142–152.
- Aurélien Beau, Denis Bossard and Sarah Gebeile-Chauty, Magnetic resonance imaging artefacts and fixed orthodontic attachments, *European Journal of Orthodontics*, 2015, 105–110.
- Blankenstein F1, Truong BT, Thomas A, Thieme N, Zachriat C., Predictability of magnetic susceptibility artifacts from metallic orthodontic appliances in magnetic resonance imaging. *J Orofac Orthop.* 2015 Jan;76(1):14-29.
- Blankenstein FH1, Truong B, Thomas A, Schröder RJ, Naumann M., Signal loss in magnetic resonance imaging caused by intraoral anchored dental magnetic materials, *Rofo.* 2006 Aug;178(8):787-93. [Article in German].
- Costa A., Appenzeller S., Yasuda C., Pereira F., Zanardi V., Cendes F. Artifacts in brain magnetic resonance imaging due to metallic dental objects. *Med Oral Patol Oral Cir Bucal* 14(6):E278-82 (2009).
- Destine D., Mizutani H., Igarashi Y. Metallic Artifacts in MRI Caused by Dental Alloys and Magnetic Keeper. *J Jpn Prosthodont Soc* 52 : 205-210, 2008.
- Eggers G., Rieker M., Kress B., Fiebach J., Dickhaus H., Hassfeld S. Artefacts in magnetic resonance imaging caused by dental material. *MAGMA* (2005) 18: 103–111.
- Elison JM1, Leggitt VL, Thomson M, Oyoyo U, Wycliffe ND, Influence of common orthodontic appliances on the diagnostic quality of cranial magnetic resonance images. *Am J Orthod Dentofacial Orthop*, 2009 Jan;135(1):8.
- Gegauff A., Laurell K., Thavendrarajah A., et al. A potential MRI hazard: forces on dental magnet keepers. *J Oral Rehabil* 1990; 17:403-410.
- Go K., Kamman R., Mooyaart E. Interaction of metallic neurosurgical implants with magnetic resonance imaging at 1.5 Tesla as a cause of image distortion and of hazardous movement of the implant. *Clinical Neurology and Neurosurgery*, Volume 91, Issue 2, 1989, Pages 109-115.
- Haacke E. M., Brown R. W., Thompson M. L., Venkatesan R. *Magnetic Resonance Imaging: Physical Principles and Sequence Design.* John Wiley, 1999.
- Hornak J. P. *The Basics of MRI. Interactive Learning Software 2008.* Edizione italiana a cura di: Larobina M. e Alfano B. - CNR - Istituto di Biostrutture e Bioimmagini, Napoli, III Edizione - Ottobre 2010.
- Jin LY1, Lin J., Influence of dental metallic materials on MR imaging, *Zhejiang Da Xue Xue Bao Yi Xue Ban.* 2009 May;38(3):328-32. [Article in Chinese].
- Lissac M., Metrop D., Brugigrad, et al. Dental materials and magnetic resonance imaging. *Invest Radiol* 1991;26:40-45.
- Lomurno G., Lucarini L., Isaza Penco S., Regi L. Gaudio S. Bracket ortodontici in corso di risonanza magnetica nucleare. *Mondo Ortodontico* 2003; 6: 445-452.
- Mattson J, Simon M. *The Pioneers of NMR and Magnetic Resonance in Medicine: The Story of MRI.* Jericho & New York: Bar-Ilan University Press, 1996.
- New P., Rosen B., Brady T., Buonanno F., Kisler J. Burt C., Hinshaw W., Newhouse J., Pohost G., Taveras J. Potenzial Hazards and Artifacts of Ferromagnetic and Nonferromagnetic Surgical and Dental Materials and Devices in Nuclear Magnetic Resonance Imaging. *Radiology* 1983; 147: 139-148.
- Nicolini C., Rigo A. *Biofisica e Tecnologie Biomediche.* Zanichelli, 1994.
- Patel A., Bhavra G. S., O'Neill J. R. S. MRI scanning and orthodontics. *Journal of Orthodontics*, Vol. 33, 2006, 246–249.
- Roth C. *MR Safety*, OutSource Inc, 1996.
- Sadowsky P., Bernreuter W., Lakshminarayanan A., Kennedy P. Orthodontic Appliances and Magnetic Resonance Imaging of the Brain and Temporomandibular Joint. *The Angle Orthodontist*, January 1988: 9-20.
- Sanders MA1, Hoyjberg C, Chu CB, Leggitt VL, Kim JS., Common orthodontic appliances cause artifacts that degrade the diagnostic quality of CBCT images, *J Calif Dent Assoc.* 2007 Dec;35(12):850-7.
- Shafiei F., Honda E., Takahashi H., Sasaki T. Artifacts from dental casting alloys in magnetic resonance imaging. *J Dent Res* 82 (8), 2003. Pp 602-606.
- Sherlock F., Crues J. High-Field-Strength MR Imaging and Metal Biomedical Implants: An Ex Vivo Evaluation of Deflection Forces. *American Journal of Roentgenology* 1988; 151: 389-392.
- Sherlock F., Curtis J. *MR Imaging and Biomedical Implants, Materials, and Devices: An Updated Review.* *Radiology* 1991; 180: 541-550.

- Starčuk Z., Bartušek K., Hubálková H., Bachorec T., Starčuková J., Krupa P. Evaluation of MRI artifacts caused by metallic dental implants and classification of the dental materials in use. *Measurement Science Review*, Volume 6, Section 2, No. 2, 2006.
- Suh J., Jeong E., Shin K., Cho J., Na J., Kim D. Han C. Minimizing artifacts caused by metallic implants at MR imaging: experimental and clinical studies. *American Journal of Roentgenology* 1998; 171: 1207-1213.
- U.S. Food And Drug Administration. Guidance for Industry and FDA Staff. Establishing Safety and Compatibility of Passive Implants in the Magnetic Resonance (MR) Environment. August 21, 2008.
- Valli G., Coppini G. *Bioimmagini*. Patron Editore, Bologna, 2005.
- Wolbarst A. B. *Physics of Radiology*. Prentice Hall, 1993.
- Borsa, Gulmanelli, Scannicchio. *Fisica biomedica*, Pavese, 1983.
- Hinshaw D. B. Jr., Holshouser B., Engstrom H., Tjan A., Christiansen E. Catelli W. Dental material artefacts on MR images. *Radiology* 1988; 166: 777-779.
- www.wikipedia.it, October 2016.
- Zahra Dalili Kajan, Jalil Khademi, Ahmad Alizadeh, Yasamin Babaei Hemmaty, and Zahra Atrkar Roushan, A comparative study of metal artifacts from common metal orthodontic brackets in magnetic resonance imaging, *Imaging Sci Dent*. 2015 Sep; 45(3): 159–168.
- Erika Ozawa, Ei-ichi Honda, Kulthida Nunthayanon Parakonthon, Hiroko Ohmori, Kazuo Shimazaki, Tohru Kurabayashi and Takashi Ono, Influence of orthodontic appliance-derived artifacts on 3-T MRI movies, Ozawa et al. *Progress in Orthodontics* (2018) 19:7.
- Maria Francesca Sfondrini, Lorenzo Preda, Fabrizio Calliada, Lorenzo Carbone, Luca Lungarotti, Luisa Bernardinelli, Paola Gandini and Andrea Scribante, Magnetic Resonance Imaging and Its Effects on Metallic Brackets and Wires: Does It Alter the Temperature and Bonding Efficacy of Orthodontic Devices?, *Materials* 2019, 12, 3971.
- Ingrid Różyło-Kalinowska, Beata Walawska, Agnieszka Predko-Engel, Elżbieta Jurkiewicz, Andrzej Urbanik, MAGNETIC RESONANCE IMAGING IN ORTHODONTIC PATIENTS Guidelines of the Polish Orthodontic Society (PTO), the Polish Medical Radiological Society (PLTR), and the Polish Dental Association (PTS), *J Stoma* 2019; 72, 1: 1–3.
- Adrienn Dobai, Fanni Dembrowszky, Tamás Vízkelety, Péter Barsi, Fanni Juhász and Csaba Dobó-Nagy, MRI compatibility of orthodontic brackets and wires: systematic review article, Dobai et al. *BMC Oral Health* (2022) 22:298.