LEOCRIL[®]

Made in Italy

RESINA ORTODONTICA AUTOPOLIMERIZZANTE ORTHODONTIC SELF-CURING RESIN RESIN ORTHODONTIQUE AUTOPOLYMERISANTE KIEFERORTHOPÄDISCHER AUTOPOLYMERISATKUNSTSTOFF

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Plaster models: Fix the various orthodontic devices (clasps, expansion screws etc..) to the model by means of sticky wax Leone R3713-00. Before preparing acrylic, dip the model into water for some minutes. To get a well made orthodontic plate, it is necessary to work on a wet, air bubble free model.

Spray-on method: After dipping of model into water for some minutes, isolate it with Leone R6320-00. Spray-on some powder (polymer) with the flask included in the package, starting from the teeth towards the centre of palate. Immediately afterwards spray on the liquid (monomer) in a sufficient quantity to moisten the powder. Proceed in constructing the plate by adding powder and liquid stepwise. Once the construction of the plate is over, spray on a light layer of powder in order to soak up the monomer excess before dipping into pressure vessel Leone T1330-00 or hydroflask Leone T1315-00. This will reduce any residual of monomer or shrinkage of acrylic.

Mixing method: After dipping of model into water for some minutes, isolate it with Leone R6320-00. Mix a ratio of 2,5 - 3 portions of powder to 1 portion of liquid into a little bowl. Keep it closed for about 3 minutes till a Plasticine consistent mixture is obtained. Spread some silicon cream on your fingers and mould the mixture on the model with a spatula (preferably moistened with monomer to avoid that acrylic sticks together) or by hand. Modelling time takes about 5 minutes.

Functional appliances: To construct functional appliances separately mould upper and lower part and joint them before dipping into Leone T1330-00 pressure vessel or Leone T1315-00 hydroflask.

Air pressure curing: It should be executed at a temperature of about 40° and under a pressure of 2,0-2,2 bar (about 2,5 Atm) for about 20 minutes. Before polymerization time is over, it is possible to interrupt the process and insert other plates or devices into the pressure vessel unless that acrylic of the other plates have remained under pressure for at least 5 minutes. Each plate shall be polymerized for at least 20 minutes in total.

A material safety data sheet in Italian and English language is available.

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