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EXPANDER/EXPANSION SCREW FOR TADs - INSTRUCTIONS FOR USE

Warning. The herewith furnished information is NOT intended for the patient, but exclusively for the medical professional, who will have to evaluate the contents in reference to every specific case to be treated, in order to furnish suitable information to the patient about the benefits and eventual side effects or complications. The expander/expansion screw for TAD must only be used by Dentists or Dental Surgeons. The hereunder instructions for use are finalized to underline the risks, the contraindications and the collateral effects of general character (a reference to all the possible circumstances of use is not possible) and well known at the moment of the placing on the market of the device. The duty of monitoring the procedures of use of the product is not ascribed to Leone; such procedures must integrally be reverted under the responsibility of the clinician. Leone takes no responsibility for the application of the device nor for possible surgical misuse of the device even neither for failure, adverse reactions or damages occurred to the patient or to the dentist following the application of the product. In any case, the use of components, accessories and instruments indicated below in accordance with the instructions is an essential condition. The expander/expansion screw for TAD performs its function through anchoring to bone or bone-dental tissue, for this reason joint use is strictly required and only with the connection/anchoring/activation devices indicated by Leone. The activation of RPE for TADs expander/expansion screw must always be carried out by the Dental Surgeon or by the Dentist (hereinafter called doctor) who decides the therapy to follow and takes care of the application of the device. **U.S.A. warning.** Caution: federal law restricts this device to sale by or on the order of a physician. **Description of the device.** The expander/expansion screw for TAD is made of stainless steel. It is a screw for the palatal expansion, included in the packages there are four TAD holders to anchor the expander on the mini-implants for palatal anchorage made up of two parts (one endosseous and a cover screw) The expander/expansion screw for TAD features two arms modelled in the dental lab and anchored to the teeth by means of orthodontic bands according to the prescription of the doctor. The activation mechanism remains completely extra mucosal. Intended Use. The expander/expansion screw for TAD is a medical device for professional use only for the palatal expansion, characterized by the anchoring to the bone by means of the mini-implants for palatal anchorage made up of two parts (one endosseous and a cover screw) and by the anchoring to the teeth by means of orthodontic bands. It must be activated with the provided Leone activation key. The expander/expansion screw for TAD is used temporarily and must be removed after orthodontic treatment has been completed. The device is single use only. It must be disposed of in the most appropriate manner and according to the laws in force. **Contraindications.** Since implant-based surgery is required, all the pathological conditions, both local and systemic, which contraindicate an intervention of oral surgery, have to be considered contraindications valid for the application of the expander/expansion screw for TAD. As information only and provided that a specific clinical evaluation of the patient is strictly recommended, the expander/expansion screw for TAD should not be used in patients with recent myocardial infarct, haematological disorders, metabolic bone disorders, neoplastic pathologies diagnosed during the last 5 years, immunosuppressed patients, previous bone radiotherapy or chemotherapy, patients with cardiovascular or metabolic problems which are not compensated, dismorphicophic patients, severe mental illness, xerostomia, diabetes or hyperglycaemia, transplanted patients, renal insufficiency, fibrous dysplasia, pharmacological treatments with anticoagulants, anticonvulsants, immunosuppressants, steroids, acute or chronic infectious illnesses, pregnancy, nursing, scarce oral hygiene, parodontitis, candidosis and pathologies of the oral mucosa, occlusion and/or articulation defects, para-functional habits, poor patient's compliance and motivation, abuse of drugs, alcohol and tobacco, acute, sub-acute or chronic maxillary osteitis. It is understood that the final choice to proceed or not with the treatment is an exclusive competence of the sanitary clinician. **Side effects.** Since implant-based surgery is required, the surgical application can be followed by edema, pain, haematomas, inflammation, limitation or alteration of the sensibility and oral functions, hypersensitivity to metal, allergic reactions. Generally, such effects are temporary and the patient should be informed about these side effects and the necessity to give immediate communication to the physician when they would occur. **Packaging and sterility.** The expander/expansion screw for TAD is non-sterile, so it must be sterilized before use. It must be sterilized according to a validated procedure which is constantly monitored and follows ISO 17665-1 "Sterilization of health care products - moist heat - Part 1: requirements for the development, validation and routine control of a sterilization process for medical devices". It is recommended to sterilize the expander/expansion screw for TAD with steam autoclave using the following parameters: temperature: 121°C (250°F), pressure: 1 atm, minimum exposure time: 20 minutes. Do not exceed 121°C (250°F), to avoid compromising functions properties of product's materials. The sterilization packaging must be appropriate for the device and conforming to ISO 11607-1 "Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems" and EN 868-5 "Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods". Never use non-sterile or pre-used the expander/expansion screw for TAD. Product must be carefully cleaned prior to sterilization. Identification and traceability. The expander/expansion screw for TAD is identified by means of a code and a lot number. This information is needed to ensure traceability of the individual product. **Final precautions.** The patient should be informed about the precautions for the period after installation of the expander/expansion screw for TAD, during and at the completion of treatment, in order to prevent complications and variations in the efficiency of the device. For instance, a good level of oral hygiene should be achieved and maintained, periodical check-ups should be performed. **Activation Procedure.** The activation of the expander/expansion screw for TAD must always be affected by the doctor who decides the therapy to follow and took care of the application of the prosthesis. Considering that the activation of the expander must be performed in a very short time, upon indications of the doctor even twice a day, he will decide whether to have the patient or patient's parent activate the prosthesis according to the home activation procedure (see the leaflet INFO82 "Informativa to activate the rapid palatal expander") for the use of the provided activation key A0556-01. The activation key A0556-01 is able to activate the device if the threshold of 50 Ncm is exceeded the key bends.

Notice for the user and/or patient: in the event of a serious accident occurring in connection with the device, immediately inform the manufacturer and the competent authority of your Country. Contact points of the European Competent Authorities: https://ec.europa.eu/health/md_sector/contact_en

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