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TADS TEMPORARY ANCHORAGE DEVICES - INSTRUCTIONS FOR USE

Warning. The herewith furnished information is NOT intended for the patient, but exclusively for the medical professional, who will have to evaluate the contents in reference to every specific case to be treated, in order to furnish suitable information to the patient about the benefits and eventual side effects or complications. Leone TADS must only be used by Dentists or Dental Surgeons. The hereunder instructions for use are finalized to underline the risks, the contraindications and the collateral effects of general character (a reference to all the possible circumstances of use is not possible) and well known at the moment of the placing on the market of the device. For the application of the product, follow the example drawings in the catalogue. The duty of monitoring the procedures of use of the products is not ascribed to Leone; such procedures must integrally be reverted under the responsibility of the clinician. Leone takes no responsibility for the application of the device nor for possible surgical misuse of the device even neither for failure, adverse reactions or damages occurred to the patient or to the dentist following the application of the product. In any case the use of exclusively original components, accessories and instruments in conformity to the instructions is strongly recommended. Since the implant stability is essential to the success of orthodontic treatment, the TADS should be stable prior to loading and implant stability should be monitored throughout the course of treatment. Depending on the host's bone anatomy, it is always indicated to place the mini-screws in the best quality bone. If inserted with a surgical motor, do not apply torques greater than 40 N-cm when inserting the mini-screw. TADS are self-drilling, however, in case of high bone density, pre-drill with one of the drills available in the orthodontic product catalogue (REF 090-1134-00, 090-1334-00). Do not place the mini-screw in contact with the roots of adjacent teeth. **U.S.A. warning.** Caution: federal law restricts this device to sale by or on the order of a physician. **Description of the product.** Leone TADS are medical devices made of titanium for surgical use (ISO 5832-3). TADS are self-tapping orthodontic mini-screws with a right-hand thread, of different lengths, supplied in sterile packaging. The fastening components (REF: 023-1702-20, 033-3203-10, 033-3203-20), which ensure the anchorage of the orthodontic device to the TADS, are medical devices made of titanium for surgical use (ISO 5832-3) and stainless steel for surgical use (ISO 5832-1), supplied in non-sterile packaging. These products are designed to be used with Leone surgical accessories and instruments, as specified in the dedicated section of the orthodontic product catalogue (REF: 080-1008-20, 083-1005-01, 083-1006-01, 156-1019-00, 083-1007-01, 080-1000-24, 080-1000-25, 080-1000-00, 080-1000-03, 043-0001-00, 043-0000-35, 083-0004-00, 053-0001-00, 083-0001-00, 083-0002-00, 083-0003-00, 070-8000-01, 080-1001-00, 080-1001-01, 080-1002-00, 080-1002-01, 080-1003-00, 156-1001-01, 156-1014-00, 080-1000-22, 080-1000-23). **Intended use.** TADS are invasive medical devices used in the oral cavity. They are devices intended to be partially implanted in the bone tissue of the oral cavity by means of a surgical procedure and removed at the end of the treatment, with the function of anchorage for orthodontic tractions and/or expanders in patients with adequate bone availability. The device is single use. **Contraindications.** All the pathological conditions, both local and systemic, which contraindicate an intervention of oral surgery, have to be considered contraindications valid for the application of the TADS as well. As information only and provided that a specific clinical evaluation of the patient is strictly recommended, Leone TADS should not be used in patients with recent myocardial infarct, haematological disorders, metabolic bone disorders, neoplastic pathologies diagnosed during the last 5 years, immunosuppressed patients, previous bone radiotherapy or chemotherapy, patients with cardiovascular or metabolic problems which are not compensated, dismorphismophobic patients, severe mental illness, xerostomia, diabetes or hyperglycaemia, transplanted patients, renal insufficiency, fibrous dysplasia, pharmacological treatments with anticoagulants, anticonvulsants, immunosuppressants, steroids, acute or chronic infectious illnesses, pregnancy, nursing, scarce oral hygiene, parodontitis, candidosis and pathologies of the oral mucosa, occlusion and/or articulation defects, para-functional habits, poor patient's compliance and motivation, abuse of drugs, alcohol and tobacco, acute, sub-acute or chronic maxillary osteitis. It is understood that the final choice to proceed or not with the treatment is an exclusive competence of the sanitary clinician. Please note that the Leone TADS and fastening components are disposable; reuse could result in cross-infection and improper functionality due to probable mechanical damage during removal. **Side effects.** Apart from the features of the utilized implant, the surgical application can be followed by edema, pain, hematomas, inflammation, limitation or alteration of the sensibility and oral functions, hypersensitivity to metal, allergic reactions. Generally, such effects are temporary and the patient should be informed about these side effects and the necessity to give immediate communication to the physician when they would occur. **Complications and interactions.** Referable complications related to Leone TADS are not known neither predictable. The surgical procedure, particularly if defectively performed, may cause chronic pain or discomfort due to the presence of the mini-screw, hemorrhagic diseases, permanent paresthesias and/or dyesthesias, bone resorption, infectious processes, damages to adjacent teeth, mini-screw migration, bending, fracture or loosening, bone fractures, mucous ulcers, gingival hyperplasia, oronasal or oronasal fistulae, aesthetical problems, nervous lesions, inhalation or ingestion of parts or components of the orthodontic implant system. In case of complications and interactions, the patient should be informed about the necessity to give immediate communication to the physician when they would occur. **Packaging and sterility.** Leone TADS are supplied gamma-sterilised in a package with a sterility indicator (the red colour of the label indicates the sterility of the product). The product is sterile if the packaging is intact (the sterility barrier consists of the tube plus the implant holder cap). The product remains sterile if packaging is intact. When stored correctly, the sterility of the product is guaranteed up to the expiry date. Do not use the product after the expiry date printed on the label. Store the product in the original packaging, in a dry place, at room temperature and away from sun light. Leone takes no responsibility in case of re-sterilized mini-screw, regardless of the method of sterilization. Never use non-sterile or pre-used TADS. To open the package and handle sterile products, please refer to the explanatory illustrations in the catalog. The sterile package must be opened in a sterile field. The accessories must be sterilised, in accordance with ISO 17665-1 "Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices", according to the following validated procedure: sterilise the accessories for orthodontic mini-screws in an autoclave using the following parameters: temperature of 121°C (250°F), pressure of 1 atm, minimum exposure time of 20 minutes. Do not exceed 121°C (250°F), to avoid compromising the properties of the material functions of the device. Use a sterilisation bag that complies with ISO 11607-1 "Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems" and EN 868-5 "Packaging for terminally sterilised medical devices - Part 5: Heat-sealable bags and tubes made on one side of porous material and on one side of plastic film - Requirements and test methods". Do not use the orthodontic device if it is not sterilised or has been used before. The product must be thoroughly cleaned before sterilisation. **Identification and traceability.** Leone TADS and fixing accessories are identified by means of a code and a lot number. This information is needed to ensure traceability of the individual product. The package is provided with an adhesive label showing these data and to be applied on the patient's record sheet. **Disposal.** Depending on the composition of the materials of Leone orthodontic and dental products, waste, generated after or before end use, requires simple or particular management to minimise its effect on the environment. Below are some of the typical methods of disposal of orthodontic and dental products: normal commercial waste, usually disposed of in landfill or, where possible, recyclable or composted, such as plastic, paper and glass packaging waste, and plastic or metal waste from laboratory products; controlled waste, usually disposed of in landfill or, where possible, which can be recycled or composted, at the end of treatment after a decontamination procedure, such as single-use devices or worn-out reusable devices; special waste, usually incinerated through specialised waste management or disposed of in landfill for hazardous waste after specialised management and pre-treated, or partially disposed of in landfill after a specialised recycling procedure, such as devices with integrated human biological residues or infected clinical waste, or expired hazardous devices that must be discarded according to the relevant materials safety data sheet. **Final precautions.** The patient should be informed about the precautions for the period after installation of the device, during and at the completion of treatment, in order to prevent complications and variations in the efficiency of the device. For instance, a good level of oral hygiene should be achieved and maintained, periodical check-ups should be performed. **For customers in the EU:** To request a copy of the Summary of Safety and Clinical Performance (SSCP), visit our website at www.leone.it. **Notice for the user and/or patient:** in the event of a serious accident occurring in connection with the device, immediately inform the manufacturer and the competent authority of your Country. Contact points of the European Competent Authorities: https://ec.europa.eu/health/md_sector/contact_en

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