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## REUSABLE SURGICAL INSTRUMENTS FOR DENTAL USE - INSTRUCTIONS FOR USE

**Warnings and guarantee restrictions.** This information is NOT addressed to the patient, but to the medical professional only, who will have to evaluate the contents with reference to each specific case, to provide adequate information to the patient about the benefits, contraindications, or complications. Reusable surgical instruments must only be used by Dentists and Dental Surgeons experienced in implant techniques as well as in the related diagnostic methods, pre-operative planning, and manufacture of prosthesis on implants. The hereunder instructions for use are finalized to underline the risks, the contra-indications and the collateral effects of a general nature (a reference to all the possible circumstances of use is not possible) and well known at the time the device was placed on the market. Nevertheless, the Leone Implant System requires the knowledge of the surgical and prosthetic procedures illustrated in the Leone Guidelines: *XCN® Implant System Surgical procedure* available on Leone Implant Catalogue or on the web site [www.leone.it](http://www.leone.it), *TADs PA and TADs VL - Insertion procedure* available in the Leone orthodontic catalogue and on the website [www.leone.it](http://www.leone.it). Considering where errors and inaccuracies in patient evaluation, pre-operative diagnosis, in the planning and in the carrying out of treatment may lead to implant failure, it is strongly recommended to attend specific training courses to with the purpose to reach a high degree of propaedeutic knowledge and practice before use. It is not ascribed to Leone the duty of monitoring the procedures of use of the product for the implant prosthetic treatment; such procedures must entirely be reverted under the responsibility of the clinician. Leone therefore takes no responsibility for the application of the device and its workmanship nor for possible prosthetic misuse of the device even neither for failure, adverse reaction or damage occurred to the patient or to the dentist following the application of the product. In any case the use of original components, accessories, and instruments only in accordance with the instructions, is strictly recommended. **Product description.** Reusable surgical instruments are medical devices made of stainless steel or titanium for surgical use. They are available in various shapes and sizes, depending on their intended use and to meet the individual anatomical variability. **Intended use.** Reusable surgical instruments are hand operated invasive medical devices intended for the implant site preparation and/or for ancillary surgical procedures. They can only be reused following a dedicated reconditioning and sterilization procedure. **Use.** Detailed information on the pre-operative planning methods, the choice of the type of instrument, the preparation of the implant site, the reopening of the implant and the precautions to be taken are contained in the specific Leone Manuals indicated above. To use the instruments maximum 20 times. **Contraindications. Side effects. Complications and after-effects.** All pathological conditions, both local and systemic, that contraindicate an intervention of oral surgery must be considered valid contraindications also for the use of reusable surgical instruments and for the application of implants and TADs. In the evaluation of the patient and in the use of the product, Dentists and Medical Surgeons must consider the contraindications, side effects, complications and aftereffects generally valid for dental surgery. Using the instrument beyond the maximum allowable number of cycles may cause damage to the device, possibly resulting in injury to the patient. **Cleaning, maintenance, and sterilization.** Leone Spa has validated two different methods for the cleaning and disinfection of reusable instruments, one manual and one automated with the use of a thermal washer-disinfector (washing and disinfecting device). If available, we recommend using the automated method. The following workflow shows the necessary activities for both validated methods. The following pages describe each activity in detail. **1. Collection and transport:** Place used instruments in a container separately from unused instruments. Prevent organic residues (blood, bone chips etc.) from drying on the instruments. To do this, clean highly contaminated instruments e.g. with gauze soaked in sterile distilled water and immerse them in sterile distilled water. Proceed with cleaning and disinfection immediately after surgery. Collect all instruments in a rigid, airtight container for safe transport to the reprocessing area. Place damaged or single-use instruments (e.g. implant carrier) in an appropriate container for separate processing before disposal. **2. Disassembly:** Disassemble instruments used in conjunction with other instruments before cleaning and disinfection. Dismantle separable instruments before cleaning and disinfecting. **3a. Manual cleaning, disinfection and drying:** Sort instruments according to material (stainless steel, titanium, aluminium, plastics) and process them separately. Follow the instructions at the end of this section for cleaning and disinfecting the ratchet and the abutment seater. Immerse the instruments in lukewarm water (<40°C/104°F) to remove organic residues. Do not use fixation agents or hot water (>40°C/104°F). Immerse the instruments in an enzymatic cleaning solution<sup>1</sup> prepared with lukewarm water. For concentration, temperature and exposure time follow the detergent manufacturer's instructions for use. Clean each instrument with appropriate brushes and pipe cleaners with soft plastic bristles under cold running water, paying attention to critical areas of the instrument such as cavities and knurls. Immerse the instruments in an ultrasonic bath with a cleaning solution.<sup>2</sup> For concentration, temperature and exposure time, follow the detergent manufacturer's instructions for use. Instruments with cutting edges must not come into contact with each other during ultrasonic treatment. Rinse instruments thoroughly with demineralised water to remove all detergent residues. Immerse the instruments in a disinfection solution.<sup>3</sup> The instruments must be completely covered by the solution. For concentration, temperature and exposure time follow the disinfectant manufacturer's instructions for use. Rinse each instrument thoroughly with demineralised water for at least 1 minute or immerse the instruments for at least 1 minute in a large quantity of demineralised water. Repeat the rinsing operation at least 3 times to remove all residues of the disinfectant. Use clean water for each rinse. Dry the instruments with compressed medical air and disposable medical wipes that do not release fibres. If necessary, repeat the cleaning and disinfection process. Ratchet and abutment seater, instruments that must not be dismantled: wash the ratchet head and abutment seater under running water to remove any residue. Thoroughly clean the outer surface of the devices with a volatile, non-viscous disinfectant. Dry the devices carefully using compressed medical air. **3b. Automated cleaning, disinfection and drying:** Sort instruments according to material (stainless steel, titanium, aluminium, plastics) and process them separately. Do not treat aluminium products, the ratchet and the abutment seater with automated cleaning and disinfection procedures. Immerse the instruments in lukewarm water (<40°C/104°F) to remove organic residues. Immerse the instruments in an enzymatic cleaning solution<sup>1</sup> prepared with lukewarm water. For concentration, temperature and exposure time follow the detergent manufacturer's instructions for use. Clean each instrument with appropriate brushes and pipe cleaners with soft plastic bristles under cold running water, paying attention to critical areas of the instrument such as cavities and knurls. In case of highly contaminated instruments, immerse the instruments in an ultrasonic bath with a cleaning solution.<sup>2</sup> For concentration, temperature and exposure time, follow the detergent manufacturer's instructions for use. Instruments with cutting edges must not come into contact with each other during ultrasonic treatment. Rinse instruments thoroughly with cold running water to remove all detergent residues. Place the instruments in a tray suitable for automated washing, close it and load it into the thermal washer-disinfector. Place the completely disassembled surgical kit and organizer in the thermal washer-disinfector. The instruments must not come into contact with each other during automated washing. Start the appropriate washing and disinfection program for the instruments. Determine the individual parameters, follow the thermal washer-disinfector manufacturer's instructions.<sup>8</sup> Start a drying cycle. Remove instruments from the thermal washer-disinfector at the end of the program. If necessary, finish drying the instruments manually with medical compressed air and disposable medical wipes that do not release fibres. **4. Inspection and assembly:** Visually inspect all instruments and accessories. Where applicable, check if the cutting edges are damaged, signs of corrosion are visible, marking is visible, the DLC coating is intact, the ISO connection is intact, the instrument is intact. Subject the instruments to a functional test; check functioning of the ratchet and the instrument coupling devices. Discard all instruments with signs of damage or loss of functional characteristics. **5. Packaging:** Check that the instruments are completely dry before inserting them in sterilization pouches. Place one or more instruments in each pouch. For instruments with cutting edges, we recommend single pouching. Close the pouch with a medical heat sealer. For steam sterilization use a sterilization pouch suitable for the device and in accordance with ISO 11607 and EN 868-5. Label each pouch with the expiry date of sterility, name of the person carrying out packaging and sterilization, name and article code of the product (if the device is not clearly visible). **6. Sterilization:** Check that the autoclave is class B, in accordance with EN 13060 and EN 285. Place the packed devices in the autoclave. Follow the autoclave manufacturer's instructions regarding placement of devices within the autoclave, paying particular attention to not overloading the autoclave and avoiding contact of the sterilization pouches with the interior walls of the autoclave and the filter area. Sterilize surgical instruments in an autoclave using the following parameters: temperature of 121°C (250°F), pressure of 1 atm and minimum exposure time of 20 minutes. Remove the packed devices from the autoclave at the end of the program. Place the packed devices on a dry, clean surface protected from air movement and cover with a clean, dry cloth for at least 10 minutes (the verification of the effectiveness of the sterilisation procedure by the user in accordance with the procedure defined above was carried out after 13 days). **7. Storage:** Store the instruments in the sterile package in a dry and closed place, such as a cart or drawer, and avoid all actions that may damage them: shocks, dropping the packages, rubbing. Follow pouch manufacturer's instructions regarding storage conditions and expiry date of sterilized products. **Notice for the user and/or patient:** in the event of a serious accident occurring in connection with the device, immediately inform the manufacturer and the competent authority of your Country.

Contact points of the European Competent Authorities: [https://ec.europa.eu/health/md\\_sector/contact\\_en](https://ec.europa.eu/health/md_sector/contact_en)

<sup>1</sup>Detergent used for validation by Leone Spa: Cidezime®

<sup>2</sup>Detergent for ultrasonic bath used for validation by Leone Spa: Cidezime®

<sup>3</sup>Disinfectant used for validation by Leone Spa: Cidex® OPA

<sup>8</sup>The following equipment, parameters and materials were used in the Leone Spa validation procedure:

Thermal washer-disinfector: Miele G 7835 CD Program: DIS. VAR. TD Detergent: Dr. Weigert ProCare Dent 10 MA

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