



The Leone expansion screws for removable orthodontic appliances derive from 80 years of extensive research, design and millions of clinical cases. The entire production process takes place in Italy with computer aided automated machines.

The screws are made of biomedical stainless steel of the highest quality and male screw threads with centesimal tolerances and total absence of surface roughness.

The bodies and the guides are made with tolerances of +/-0.01 mm which ensure stability and precision. The wide range available allows the production of any type of removable device.


excellent  
for stability  
and absolute  
precision



■ built-in screw stop at the stroke end

■ safe retention into acrylic

■ anti-slippage screw



**REMOVABLE**

I.P.

Italian  
style

IN ORTHODONTICS  
SINCE 1934



LeoneAmerica  
1250 Stellar Drive | OXNARD, CA 93033-2404  
805-487-9860 | 800-235-4639 | info@leoneamerica.com | www.leoneamerica.com  
LeoneAmerica sole distributor of Leone S.p.A. orthodontic products in USA, Canada and Puerto Rico

IN-09-19



Orthodontics and Implantology

Italian  
style

ORTHODONTIC  
SCREWS



Leone palatal expanders have been known worldwide for more than 50 years and are now the most complete range of expanders in the market.

The 620 series is a superior product made in various dimensions with many advanced features including laser welding of the extension arms into the body of the expander for greater strength and stability.

The automated manufacturing process and quality control with the highest standards allow us to confidently say that we are offering the best palatal expanders in the world today.


A0620  
the best-selling  
expander  
ever in the world



laser-etched ■

made entirely of  
biomedical  
stainless steel ■

laser welding  
of the extension arms  
into the screw's body ■



**FIXED**

**RAPID EXPANDER**

A0620-08  
8 mm

A0620-09  
9 mm

A0620-11  
11 mm

A0620-13  
13 mm

A0620-15  
15 mm

**ZOOMRAX RAPID EXPANDER**  
WITH TELESCOPIC GUIDES

A2620-07  
7 mm

A2620-10  
10 mm

A2620-12  
12 mm

**EAGLE RAPID EXPANDER**  
WITH ORTHOGONAL ARMS

A0630-08  
8 mm

A0630-10  
10 mm

A0630-12  
12 mm

**stealth MICRO EXPANDER**  
FOR PALATAL SUTURE

A0621-08  
8 mm

A0621-11  
11 mm

**stealth SINGLE ARM MICRO EXPANDER**  
EASY ACCESS

A1621-08  
8 mm

A1621-11  
11 mm

A1621-13  
13 mm

**EXSPIDER™ RAGNO FAN-TYPE**  
RAPID EXPANDER

A0625-09  
9 mm

A0625-07  
7 mm

**TORKO™ EXPANDER**  
FOR LOWER ARCH

A0623-08  
8 mm

A0623-11  
11 mm

**LEAF EXPANDER® WITH MEMORIA®**  
SPRING

A2703-06 2 springs 450 g approx.  
6 mm

A2704-06 2 springs 900 g approx.  
6 mm

A2703-09 3 springs 450 g approx.  
9 mm

A2704-09 3 springs 900 g approx.  
9 mm

**LEAF self EXPANDER® WITH MEMORIA®**  
SPRING

A2705-06 4 springs 450 g approx.  
6 mm

A2706-06 4 springs 900 g approx.  
6 mm

A2705-09 6 springs 450 g approx.  
9 mm

A2706-09 6 springs 900 g approx.  
9 mm

**MULTIFUNCTIONAL**  
ORTHODONTIC SCREWS

A3626-11  
11 mm

A3628-11 DX  
11 mm

A3629-11 SX  
11 mm

**POP SCREW**  
FOR UPPER ARCH

A4805-14G  
A4805-14E  
A4805-14R  
A4805-14V  
A4805-14T  
A4805-14A

**POP SCREW**  
UNIVERSAL

A4800-12G  
A4800-12E  
A4800-12R  
A4800-12V  
A4800-12A  
A4800-12A

**STANDARD SCREW**  
FOR UPPER APPLIANCES

A0805-11  
7 mm

A0805-14  
9 mm

**MEDIUM SCREW**

A0803-11  
7 mm

A0803-14  
9 mm

**UNIVERSAL SCREW**

A0800-10  
6 mm

A0800-12  
8 mm

**SUPER MICRO SCREW**

A0890-11  
6 mm

**MICRO SECTIONAL SCREW**

A0801-00  
6 mm

A0801-10  
4 mm

**MICRO SECTIONAL SCREW**  
WITH STRAIGHT U-SHAPED GUIDE PIN

A0801-12  
4 mm

A0801-13  
6 mm

A0801-14  
4 mm

**MICRO SECTIONAL SCREW**  
WITH BENT U-SHAPED GUIDE PIN

A0802-00  
4 mm

A0802-13  
4 mm

**SUPER MICRO SECTIONAL SCREW**

A0891-09  
3 mm

A0891-10  
3 mm

**THREE DIRECTIONAL SCREW ANATOMICAL**

A0930-14  
3 mm

A0930-16  
4 mm

3 activation points

**THREE DIRECTIONAL SCREW ANATOMICAL**

A0931-14  
3 mm

A0931-16  
4 mm

2 activation points

**THREE DIRECTIONAL SCREW ANATOMICAL FOR LOWER ARCH**

A0932-13  
2.5 mm

**FAN-TYPE SCREW**

A0515-14  
9 mm

A0516-14  
9 mm