

IMPLANT SYSTEM

XCN[®]

DEMONSTRATION AND INFORMATIVE MATERIAL



SURGICAL DEMONSTRATION KIT

- to simulate the most relevant surgical and prosthetic phases of the XCN® implant system
- intended for demonstration use only

Pack content:

- 1 hemi-mandible
- 1 Classix implant Ø 4,1 mm (non-sterile) 10 mm long with cover cap and 5 cover cap refills connection 3,0
- 1 instrument for cover caps
- 1 Standard healing cap connection 3,0 GH 3 mm
- 1 Standard transfer connection 3,0
- 1 surgical hand screwdriver
- 1 hex head extractor
- 1 straight Standard Basic abutment connection 3,0



REF **106-0002-00**

HEMI-MANDIBLE

- made of polyurethane
- presents a hole for the placement of a Classix implant Ø 4,1 mm, 10 mm long

Pack of 1



REF **106-0001-00**

DEMONSTRATION JUMBO DENTAL IMPLANT

- made of aluminium
- scale 5:1 reproduction of straight Standard Basic abutment connection 3,0 and Classix implant Ø 4,1, 10 mm long
- use the rod included in the package for removal of the abutment from the Jumbo implant. A hole is present in the back of the implant for the rod

Pack of 1



REF **106-0003-00**

IMPLANT CARD

- supplied with each implant
- for unique identification of the placed implant
- for delivery to the patient
- contains essential information in case of dental care away from home or abroad



PATIENT INFORMATION BROCHURE

- available on request

Pack content: 50 brochures in a carton box

